

Life Pen

Your portable magnetic field generator



Life Pen is a magnetic field generator device that acts on energy imbalances that generate physical problems of sciatica, muscle and joint pain. The principle of operation of the Life Pen is very similar to that of acupuncture, with the considerable difference that you can use it directly and you do not have to touch your body with needles.

An alternative pain relief

The Life Pen is excellent for the treatment of acute and chronic physical pain caused by arthritis, rheumatism, osteoporosis, sciatica, migraine, sports injuries and muscle, joint and back pain. Lightweight and portable: the Life Pen is small and fits easily into your pocket or purse. Take it wherever you go and treat yourself to a session whenever you want. Use it while travelling or when you need quick pain relief.

Technology

Our technology has distant historical origins. In China, it was used many years before the birth of Jesus Christ. Some references to this can be found in the works of Plato, Homer and Aristotle. Many experts and scientists, in various countries of East and West, have studied and conducted experiments on biomagnetism, but it was only after 1940 that world research began the scientific phase that produced texts that systematically validated the beneficial effects on plants, animals and humans.

The effects of our technology are not magical or miraculous but are well anchored in science because the cells of the human body are small natural magnets, with a North and South pole. Moreover, all the functions of the organism are regulated by electrical and magnetic impulses, without which man could not live.

How to use Life Pen

The pen can be used in 3 ways:

- By directly indicating the pain area;
- Use the reflexology points on your hands or feet to hit the deepest areas of pain;
- Place in a glass or jug of water and leave for at least 30 minutes to energize the water.

At first you may get a tingling sensation on your tongue.

Many people find that they have more energy and require less food.


Used for muscle aches, joint pain, toothache, wound pain, burns, bites, stings, back pain.

Can be used together with other **SoftMAG** devices.



There are eight essential points in the hand. Here's where you can find them, what they do, and how you can manipulate them for your benefit:

Heart 7




The heart 7 point can be found in the crease of your wrist. It lies in line with the space between your ring and pinkie finger.

There's a bone directly in line next to this point. Reflexologists claim that applying Life Pen on this spot may protect against **anxiety, insomnia, heart palpitations, and depression.**


Small intestine 3

The small intestine 3 point is located at the outside part of your hand, right on the edge. The point lies in the depression of your hand just below your pinkie finger.

Applying Life Pen on this point is thought to relieve **neck pain, earaches, and headaches that occur at the back of your head**




Lung meridian



You can find your lung meridian point by looking at the edge of your hand. It runs from the tip of your thumb down along the side of your hand, ending just beneath the crease of your wrist.


Approach our Life Pen by making a circular motion from a distance of about two inches from one of the 3 spots until it feels better. This may help relieve any **cold symptoms**, such as **chills**, **sneezing**, **runny nose**, and **sore throat**.

Inner gate point



The inner gate point can be found in a special spot on the crease of your wrist. Hold out your hands as if you're receiving a gift, with your wrists together. Take one hand and feel about 3 centimeters out from where your wrists were touching in the middle. Approach Life Pen to this point and make circular motion from a distance of about two inches from the above mentioned point. This is believed to **improve your digestion** and **relieve any nausea or stomach pain**.

Outer gate point




The outer gate point can be found between two tendons on the backside of the arm. Place three fingers from your other hand above your wrist. Approach Life Pen to this point and make circular motion from a distance of about two inches from the above mentioned point to potentially **give your immune system a quick boost and a rush of energy.**


Wrist Point 1

You can find your wrist point 1 on your wrist. Run a finger down your pinkie to your wrist crease, keeping it in line with your finger. Now you've found the point.

Approach Life Pen to this point and make circular motion from a distance of about two inches from your wrist point 1. This may **help promote happiness by regulating your emotions.**




Base of the thumb point



The base of the thumb point is located on your wrist. Run a finger down your thumb, to your wrist crease at the base of your thumb. Approach Life Pen to this point and make circular motion from a distance of about two inches to **help alleviate respiratory and breathing problems.**


Hand valley point

The hand valley point can be found in the firm skin between the thumb and index finger. Approaching Life Pen to this point may **help reduce stress**, as well as **alleviate migraines, toothaches, shoulder tension, and neck pain.**



Ten dispersion

Ten dispersion points are located on the tips of each finger. Approaching Life Pen to this point could help relieve some common flu symptoms, such as a high fever or a sore throat.



Four seams

A person can find the four seams on the inside of each of the large joints in the index finger, middle finger, ring finger, and little finger.


Approach Life Pen to this point and make circular motion from a distance of about two inches to **help treat digestive problems, especially in children.**



An mian

An Mian are traditional points for treating insomnia and also to reduce anxiety, vertigo, and headaches.


The An Mian points are on either side of the neck. To find them, place a finger behind each earlobe, and move the fingers just behind the bony protrusion. Approach Life Pen to this point and make circular motion from a distance of about five inches. Using these points in combination with others may help treat insomnia caused by depression.



Point SP6

Point SP6, also called **San Yin Jiao**, may help with **insomnia**, **menstrual cramps**, **urinary issues**, and some other **pelvic problems**.

To access this point, find the highest point of the ankle on the inside of the leg. Beginning at the top of the ankle, measure four finger-widths up the leg. Apply Life Pen just behind the bone above the ankle and make circular motion from a distance of about two inches for 3 minutes,




Point LV3

Point SP6, also called Tai Chong, may help with unexplained insomnia, as well as stress and anxiety-related sleeplessness. Find it by locating the spot where the skin of the big toe and the next toe connect. Apply Life Pen for 3 minutes.



Point KD3

Point KD3, also called Taixi, can help insomnia, This point is located just above the heel on the inside of the foot. Applying Life pen into the use of KD3 and HT7 found that magnetic field therapy on these points improved sleep quality in middle-aged and older adult participants with hypertension. It also helped lower their blood pressure to healthy levels.



Yin Tang



The Yin Tang point is in the center of the eyebrows, just above the nose. Apply Life Pen and make circular motion from a distance of about two inches for 3 minutes. may help relieve insomnia and other issues, including: fear, agitation, restlessness.

Diet and supplements do play an important role too. Consider including ActiveBlast Nano-patches to your sessions.

Life Pen

Your portable magnetic field generator

How it's made

Life Pen is made of PLA which is an innovative bioplastic. It is a material derived from the transformation of sugars present in corn, red beet, sugar cane and other natural and renewable materials and not derived from oil (unlike traditional plastic).

Materials and production: 100% Italy.

Warranty: 3 years

DO NOT point to the brain. DO NOT use on or near implanted devices.



www.green-shield.eu